



Recombinant Protein mxiG (mxiG)

Product Code	CSB-CF363182SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0A221
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MSEAKNSNLAPFRLLVKLTNGVGVDEFPLYYGNNLIVLGRITETLEFGNDNFPEN IIPVTD SKSDGIIYLTISKDNICQFSDEKGEQIDINSQFNSFEYDGISFHLKNMREDKSRG HILNG MYKNHSVFFFFAVIVVLIIFSLSLKKDEVKEIAEIIDDKRYGIVNTGQCNYILAETQ ND AVWASVALNKTGFTKCRYILVSNKEINRIQQYINQRFPPINLYVLNLVSDKAELL VFLSK ERNSSKDELDKLNALIVEFPYIKNIKFNLYSDHNARGDAKGIFTKVNVQYKEI CENNK VTYSVREELTDEKLELINRLISEHKNIYGDQYIEFSVLLIDDDFKGKSYLNSKDSY VMLN DKHWFFLDKNK
Source	in vitro E.coli expression system
Target Names	mxiG
Protein Names	Recommended name: Protein mxiG
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.