



Recombinant Human Transmembrane protein 183A (TMEM183A)

Product Code	CSB-CF023767HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8IXX5
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MARGPGPLGRPRPDTVAMPKRGKRLKFRAHDACSGRVTVADYANSDPAVVR SGRVKKAVA NAVQQEVKSLCGLEASQVPAEEALSGAGEPCDIIDSSDEMDAQEESIHERTVS RKKKSKR HKEELDGAGGEEYPMDIWLLLASYIRPEDIVNFSLICKNAWTVTCTAAFWTRLY RRHYTL DASLPLRLRPESMEKLRCLRACVIRSLYHMYEPFAARISKNPAPIESTPSTLKNS KCLLF WCRKIVGNRQEPMWEFNFKFKKQSPRLKSKCTGGLQPPVQYEDVHTNPDQD CCLLQVTTL NFIFIPVMGMIFTLFTINVSTDMRHHRVRLVFQDSPVHGGRKLRSEQGVQVILD PVHSV RLFDWWHPQYPFSLRA
Source	in vitro E.coli expression system
Target Names	TMEM183A
Protein Names	Recommended name: Transmembrane protein 183A
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.