



# Recombinant Bacillus licheniformis UPF0754 membrane protein BLi01057/BL02871 (BLi01057, BL02871)

<b>Product Code</b>	CSB-CF734282BQU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q65LU9
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Bacillus licheniformis (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
<b>Sequence</b>	<p>MYVFGIFAVMIAV GALIGAVTNHFAIKMLFRPYKAIYIFGKRV PFTPGLIPKRRDE  LARQ  MGQMVTGHLLTTEGLKKRLASDAVKSQAVQVGERLLARLSQSTATVEEALES  GISNPAQ  KADRAVSRLADEKLTAFL EAYENEPLK KLFPLEAQEKLKEKIPMVSSYILSRAVS  YFESD  EGKERLGHMIDDFLKERGMLGSMVQMFLGNSSLVDRVQPEIVKFLKNGETAG  LLQDLLEN  EWDKLKEYTFKEADDKWNLKPLIFDLKEKLLKRFSLQPFFEKTIGSSISSFEQDI  ALRLP  QMADRLLLEEAGRRRLDQALKQLELEQIVKEQVDNFPVERLEEMVLSISKREFKMI  TYLGGL LGGIIGAVQAIFVILI</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	BLi01057
<b>Protein Names</b>	Recommended name: UPF0754 membrane protein BLi01057/BL02871
<b>Expression Region</b>	1-376
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.