



Recombinant Bacillus subtilis UPF0754 membrane protein yheB (yheB)

Product Code	CSB-CF522472BRJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O07543
Product Type	Transmembrane Protein
Immunogen Species	Bacillus subtilis (strain 168)
Sequence	MGIAGTFIFMIVIGAAIGAVTNHLAIQMLFRPYKAYYLF GKRV PFTPGLIPRRRDE LAKQ MGLMVVNHLLTPEGIKKRLVSDAAKTQALRVGEQLIQKLSLSEVTVKEALEKAG MKRPEK AADAWISSWTDDKLHELFRQYGDQSLKELVPIEVQE KLEEKIPMISGYILSR SV RYFESD EGKIRLGNMIDDFLKERGMLGSMVQLFLGNSSLADRVLPELLKFLRNEETNKL SDLLKN EWGKLREYTFNEADEKWNALIFSLKRRVLQAFSTAPFFNNTIGTLTVRYESE LTQQML PALLDKLLLEGISSNLESVLKRLRLEEIVKEQVDQFPVERLEEMVLSISKKEFKMIT YLGGLGGIIGAIQALFVILF
Source	in vitro E.coli expression system
Target Names	yheB
Protein Names	Recommended name: UPF0754 membrane protein yheB
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.