



Recombinant *Oryza sativa* subsp. *indica* UPF0496 protein 3 (OsI_009784)

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-CF383699OFF |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A2XCJ1 |
| Product Type | Transmembrane Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. <i>indica</i> (Rice) |
| Sequence | MGATFRCFGGCVKPDDQQVHEPKKVVPSSSDFFREEYTSAFRTESYNDFW ARVLDITLA HGAALVPRHGGGGGCAASKRLPSYRLFAEHLLEPDQRAVAAALASPRGSRLR PDVRGLLA AYYAETANASFLCSHLLKDIEHIRLYRPLKHTLRKLASDVGVSGLADVSAALG QPFTAL AASQGRLREVQAGSGDLLRGLDAGRKKARHRIRSVARLRRALSVSFVTAVAV VAVVGACI GVHILAAFAAFPMMSPAWLGERFFSGRAARRALVQLEAAAKGTYILNRDMETI SRLVARV RDEGEHMVALLRLCVEHRPAAGAGGKGRVLVQEVLRQLSKNEESFRQQLDELE EHLFLCFM TINKARIMVMNFMAAAAR |
| Source | in vitro <i>E.coli</i> expression system |
| Target Names | OsI_009784 |
| Protein Names | Recommended name: UPF0496 protein 3 |
| Expression Region | 1-378 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |