



Recombinant *Saccharomyces cerevisiae* Protein COS8 (COS8)

Product Code	CSB-CF340000SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38723
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MKENEVKDEKSVDVLSFKQLEFQKTVLPQDVFRNELTWFCYEIYKSLAFRIWM LLWLPLS VWWKLSSNWIHPLIVSLLVFLGPFVVICGLSRKRSLKQLIQFCKEITEDTP SSDPH DWEVVAANLNSYFYENKTWNTKYFFFNAMSCQKAFKTTLLEPFSLKKDESAK VKSFKDSV PYIEEALQVYAAGFDKEWKLFNTEKEESPFDLEDIQLPKEAYRFKLTWILKRIFN LRCLP LFLYYFLIVYTSGNADLISRFLFPVVMFFIMTRDFQNMIRMIVLSVKMEHKMQFLS TIINE QESGANGWDEIAKKMNRYLFEKKVWNNEEFFYDGLDCEWFFRFFFYRLLSLK KPMWFASL NVELWPYIKEAQ SARNEKPLK
Source	in vitro E.coli expression system
Target Names	COS8
Protein Names	Recommended name: Protein COS8
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.