



Recombinant Uncharacterized protein C05B5.2 (C05B5.2)

Product Code	CSB-CF341511CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34290
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MLNIVLIIGLLAIFNTSSASNDVCHVTRKPLMLPAPGDVYKKAVQFYSNITAPRS TSVLA PVMSSLEVYINTTTTSAFAPAQSIKVADILEEDADAIRVKSIRMAGFIAQCIIFLFV YTI VTMDVEIWKINMDWLKIQYFQHFEDSAAEVPVFKLYMAREIQTCP LPARQNV ILSLIKH IAFIGFQFWLTKPPVARRMTFLKLAIQVLRIPFLFFIAFRAFIIPKFIFLETGSYVAT GL FSAFHLTILWIHMNPKSLKSLYAYIALTVLMLTGWDFVFRKHCFLLKKGFLVA MGIAG NVPNSREVTMIIELDMLS LTLTLVLLYLQLKKSVAEEFVFPKQPKKSIRGPVCFT NTLTE DNEGYFENIYGDLIFTSRTSNTFSI
Source	in vitro E.coli expression system
Target Names	C05B5.2
Protein Names	Recommended name: Uncharacterized protein C05B5.2
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.