



Recombinant Arabidopsis thaliana RING-H2 finger protein ATL43 (ATL43)

Product Code	CSB-CF704855DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5EAE9
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	DNHTAVVITTSDTPPPLPPPSPPPRHNF TSSLMPGIAVVIAVLTAFFSLTFLLLLY VKHC KRRNGSVYVNHQPQRFAITRYGGGYNGGGVVGGRKNSGIDRSVIESLPVFRF GALSGHKD GLECAVCLARFEPTEVLRLLPKCKHAFHVECVDTWLDHAHSTCPLCRYRVDPE DILLIGDC NSWFELQFSKDESNSVNNNPPGLTRFIPVSRISGRHSSAGERASRLNEIRTSS SYKSNPM SFRRSLDSSLKVNDAGEEKSESVAVNCLDRLQRKDGLLLIPNRESFEGRFEHR IISGGN RDDQRWSEVRPSDLLYL RSEMILSDCRKLA AEGGRDVINGRSVSELTGIERR RRWGGEPRQRQATAVISRWLAWSHRASASSIV
Source	in vitro E.coli expression system
Target Names	ATL43
Protein Names	Recommended name: RING-H2 finger protein ATL43
Expression Region	23-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.