



Recombinant *Saccharomyces cerevisiae* Monopolar spindle protein 2 (MPS2)

| | |
|--------------------------|---|
| Product Code | CSB-CF345263SVG |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P53159 |
| Product Type | Transmembrane Protein |
| Immunogen Species | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) |
| Sequence | MSNGAFDAIFEYAWGQIDKPISGDFIYGKDLPKLIEIENIFQKAQKSGSYELRLP LFSE INKDLFRTFSNTKTFKIHKEEFDDIFFNLVNHPLREILENAFIGVDSIPSDFIVSM NLN SPSKFLVENKNKNTEGAGISTPRKKLTESPIKLLSRNNIGKALEVQVEELKRELT AKQSL LQENERQVSELKIRLETYQEKYASIQQRFSDLQKARQVEDNQNSSRTSDPGSP LVTGIDQ KAILEEFRRRLQRQTDTISFLKDQIRRERGLNCSNDKVSHSKRKHATTDGDGT FKNFISA VPSNIWVKATIRIIVCFALLAGVLPYIRKYVYAHDTPSQNSRLQLSWWENSGILS KIVWF FEDQTDLETEYRSNANVDDAYS RVFGI |
| Source | in vitro <i>E.coli</i> expression system |
| Target Names | MPS2 |
| Protein Names | Recommended name: Monopolar spindle protein 2 |
| Expression Region | 1-387 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |