



Recombinant *Saccharomyces cerevisiae* Monopolar spindle protein 2 (MPS2)

Product Code	CSB-CF509830SVN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C8Z8H2
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Sequence	MSNGAFDAIFEYAWGQIDKPISGDFIYGKDLPKLIEIENIFQKAQKSGSYELRLP LFSE INKDLFRTFSNTKTFKIHKEEFDDIFFNLVNHPLREILENAFIGVDSIPSDFIVSM NLN SPSKFLVENKSKNTEGAGISTPRKKLTESPIKLLSRNIGKALEVQVEELKRELT AKQSL LQENERQVSELKIRLETYQEKYASIQQRFSDLQKARQVEDNQNSSRTSDPGSP LVTGIDQ KAILEEFRRLQRQTDITISFLKDQIRRERGLNCSNDKVSHSKRKHATTDGDGT FKNFISA VPSNIWVKATIRIIVCFALLAGVLPYIRKYVYAHDTPSQNSRLQLSWWENSGILS KIVWF FEDQTDLETEYRSNANVDDAYS RVFGI
Source	in vitro <i>E.coli</i> expression system
Target Names	MPS2
Protein Names	Recommended name: Monopolar spindle protein 2
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.