



Recombinant *Saccharomyces cerevisiae* Monopolar spindle protein 2 (MPS2)

Product Code	CSB-CF517207SVI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	E7KCH6
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain AWRI796) (Baker's yeast)
Sequence	MSNGAFDAIFEYAWGQIDKPISGDFIYGKDLPKLIEIENIFQKAQKSGSYELRLP LFSE INKDLFRTFSNTKTFKIHKEEFDDIFFNLVNHPLREILENAFIGVDSIPSDFIVSM NLN SPSKFLVENKSKNTEGAGISTPRKKLTESPIKLLSRNNIGKALEVQVEELKRELT AKQSL LQENERQVSELKIRLETYQEKYASIQQRFSDLQKARQVEDNQNSSRTSDPGSP LVTGIDQ KAILEEFRRRLQRQTDTISFLKDQIRRERGLNCSNDKVSHSKRKHATTDGDGT FKNFISA VPSNIWVKATIRIIVCFALLAGVLPYIRKYVYAHDTPSQNSRLQLSWWENSGILS KIVWF FEDQTDLETEYRSNANVDDAYS RVFGI
Source	in vitro <i>E.coli</i> expression system
Target Names	MPS2
Protein Names	Recommended name: Monopolar spindle protein 2
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.