



Recombinant Mouse Mitochondrial carrier homolog 1 (Mtch1)

Product Code	CSB-CF745123MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q791T5
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MGASDPEVAPWAPGGAAGMAGAGAGAGARGGAPAGVEARARDPPPAHRA HPRHPRPAAQP SARRMDGGPGAPGSGDNAPTTEALFVALGAGVTALSHPLLYVKLLIQVGHEP MPPTLGTN VLGRKVLVLPSSFFTYAKYIVQVDGKIGLFRGLSPRLMSNALSTVTRGSMKKVFP PDEMEQ VSNKDDMKTSLKKVVKETSYEMMMQCVSRMLAHPLHVISMRCMVQFVGREA KYSGLSSI GKIFKEEGLLGFFVGLIPHLLGDVVFLWGCNLLAHFINAYLVDDSVSDTPGGLG NDQNPG SQFSQALAIRSYTKFVMGIAVSMLTYPFLLVGDLMVNNCGLRAGLPPYSPVF KSWIHCW KYLSVQQQLFRGSSLLFRRVSSGSCFALE
Source	in vitro E.coli expression system
Target Names	Mtch1
Protein Names	Recommended name: Mitochondrial carrier homolog 1 Alternative name(s): Mitochondrial carrier-like protein 1
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.