



Recombinant *Schizosaccharomyces pombe* Meiotically up-regulated gene 89 protein (mug89)

Product Code	CSB-CF636971SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1MTQ5
Product Type	Transmembrane Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MPEALNENVSDTASNGPVAKTRAPPNTSFRQQRIKSWQPLLTPKIVLPLFFVL GIIFGPL GGGLLYASSIVQELVVDYTDCE TLASYDEFSAVPSKKYTASFDQSGTIGFDKES TYWKLE KILDKDLMDVNYCIIRFTVPSVLKAPIFIYYRLTNFFQNHRRYAKSVDEKQLQG VALTA DEVKGGNCFPLEVNEDDKPYYP CGLIANSLFNDFSSLRLLDDNSVYTFSTKNI AWASDK RRFLKTNYSRDDVAPPPNWVLRYPDGYTESNMPDLSTMENLQVWMRTAGLP TFSKLAMRN DNDDIFPGTYEIKIGLFFPVKSFDGTKSLVLTTTRSVLGGKNPFLGIAYIVVSAVCV VLGT VFTLRHFIRPRKLADHRYLNWDSEENNLAPHLSDRP
Source	in vitro E.coli expression system
Target Names	mug89
Protein Names	Recommended name: Meiotically up-regulated gene 89 protein
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.