



Recombinant Schizosaccharomyces pombe Uncharacterized protein C1705.03c (SPAC1705.03c, SPAC1F2.01, SPAC23H4.19)

Product Code	CSB-CF515613SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O13960
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	ASNCSSGPYNISAQGTLDLNSCTVLNGDLYISDAGNSGITTLLTVNGIESVQGD VVVSDG QYLTSLSFPSLKNVSGAFNVNMMIRMNNLATPELTSVGSLLNLAVLPNLQELQF NAGLSDS DSVVIDDTQLQAIDGISLDSVTTTFQVTNNRYIQEITMEGLESAQNIQISANSKGV SVNFS KLSNVTTATFDGISNVFIGNLKSAAGNLYFSNTTLDNISVPYLTEIGQSFAVLYSP ELTS LNFPNLTTVGGGFVINDTGLTSDGFPVISEIGGGLVLLGNFSSIDMPDLSDVKG ALTVE TKATNFTCPWSNDDSVIKGDDFTCQGSVATISATSSYDLSSTVSATSGSATSA TGSATTT SYSSDSSASSSSSSSHESSAASNGFTAGALVLGSLVAALAM
Source	in vitro E.coli expression system
Target Names	ecm33
Protein Names	Recommended name: Uncharacterized protein C1705.03c
Expression Region	20-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.