



Recombinant Schizosaccharomyces pombe Uncharacterized J domain-containing protein C17A3.05c (pi041, SPBC17A3.05c)

Product Code	CSB-CF515564SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O13633
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MESSREEALRCIGLARKYLNAGEYDKALKFANKSLRIHATTEGEDFVKQIEKQK LTGSPN PQATPKQENKSNFFSEKQSVRENGNSSAGEKKQKWTSEQHLLVQKIIKYKNH QYYEILD KKTCTDTEIKKSYKKLALQLHPDKNHAPSADEAFKMVSKAFQVLSDPNLR AHDRTGMDP ESRASAASSSFSSNAGGHPGFSAYPQANMSPEDLFNSFFGDQFFSGPGTFFF GGPGIRV HQFGGRPRNFARRQQAQDMPPKSIFYQLLPLIVVILFAFLSNFSWSDSTSVNT RYSFQQN YKYTVPRTTAKHNIPYYMSQKDLKSSRDIRRLNEKVEHTYQTQNVHNA CLRE QQIKEDE IRRAQGWFFPDKEALKKAKELRLPNCEELNRLGYRTYSNSYYF
Source	in vitro E.coli expression system
Target Names	pi041
Protein Names	Recommended name: Uncharacterized J domain-containing protein C17A3.05c
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.