



Recombinant Danio rerio Tripartite motif-containing 13 (trim13)

Product Code	CSB-CF024450DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q503I2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MELLEEDLTCPICCCFLFEDPRVLPCSHSFCKKCLEGILDGNSPTWRPPFKCP TCRKETV HNGIASLQVNYSLRGIVEKYNRIRVMRMSQCRVHSGQPLNIFCATDLKLICGF CATTGD HKGHKFCALEEAYEREKLAFFELFRVVEGWKGAEVHSCLESLES AKKKALERV SRDADRV SEYFDKLLRTLHHRSEILSDLETCLKLAVMQTFDPEINRLRSAL EEQRRALNIAE SFRSL SDPLTFLQQMQDFREKLRVIQGTPLPSRTDMDVLSLSALQSFDVKEWDRVRLG QVDKLCAP YESSAYLSSLPPAAAPRFTRVMWRVVLVVCACL PALNFLPSDCLALS FQDKVV ALGGFSL PSPGEIVRWLGFCWKEAASIC TLLTELCRNCMLDLINTT SDFIS
Source	in vitro E.coli expression system
Target Names	trim13
Protein Names	Recommended name: Tripartite motif-containing 13
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.