



Recombinant Rat Protrudin (Zfyve27)

Product Code	CSB-CF744070RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6P7B7
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MQSSDRDLSGPEASPSVMPEVLSECSAPATKSTAFDLFNLVLSYKRLEIYLEPL KDAGDG VRYLLRWQMPLCSLLTCLGLNILFLTLNEGAWYSVGALIISVPALLGYLQEVCR AQLPES ELMRRKYHSVRQEDLQVRVLSRPEAVAEVKSFLIRLEAFLARLCYTCESAYRV LHWENPV VSSQFYGALLGMVCMYLLPLCWVLLALLNSTLFLGNGEFFRVVSEYRAQLQRR MSPKQEE CVCEGSALQDAGGRAVVDSTPAPTPTEDLTPGSVEEAEEAEPDEEFKDAIEE DDEGTPC PAEDELTMQDNGFLSKNEVLRKSKLTERLRKRYPTNNGNCAGCAATFSV LKKRRSCS NCGNSFCSRCCSFKVPKSSMGATAPEAQRETVFVCASCNQTLISK</p>
Source	in vitro E.coli expression system
Target Names	Zfyve27
Protein Names	Recommended name: Protrudin Alternative name(s): Zinc finger FYVE domain-containing protein 27
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.