



Recombinant Danio rerio E3 ubiquitin-protein ligase MARCH4 (41337)

Product Code	CSB-CF606683DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0P496
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	RRPMLRRQGKQKGRCCVLFSDLEVFLLRPPTPSASPPAFTPMNELNAEGNAT SSATESHS LANGHYQPVTGEEALDTRGPDDWTHSVVDPRTLDCSSSEDCSKEKLDLERL SLNSCTDS GVRTPLCRICFQGPEQGEELLSPCRCSGVRCTHEPCLIKWISERGSWSCELCY YKYQVIA ISTKNPLQWQAISLTVIEKVQIAAAVLGSLFLIASISWLWSSLSPSAKWQRQDL LFQIC YAMYGFMDLVICIALIVHEGPSVFRIFNRWQAVNQWQKVLNYDKVKDNEDHQK TGATFRTL SLPLTHRMGQSGPEGEPSTSTSSLMAAAAAAAGVTPTTNSVPPAAGATTE PQDSSEPS NGQPSLPDHHCAYNILHLLSHLRQQEPRGQTSNSNRELVMRVTTV
Source	in vitro E.coli expression system
Target Names	march4
Protein Names	Recommended name: E3 ubiquitin-protein ligase MARCH4 EC= 6.3.2.- Alternative name(s): Membrane-associated RING finger protein 4 Membrane-associated RING-CH protein IV Short name= MARCH-IV
Expression Region	17-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.