



Recombinant Enterobacteria phage f1 Attachment protein G3P (III)

Product Code	CSB-CF304787ECH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P69169
Product Type	Transmembrane Protein
Immunogen Species	Enterobacteria phage f1 (Bacteriophage f1)
Sequence	AETVESCLAKPHTENSFTNVWKDDKTLDRYANYEGCLWNATGVVVCTGDET QCYGTWVPI GLAIPENEGGGSEGGGSEGGGSEGGGKTPPEYGDTPIPGYTYINPLDGTYP GTEQNPAN PNPSLEESQPLNTFMFQNNRFRNRQGALTVYTGTVTQGTDPVKTYQYTPVS SKAMYDAY WNGKFRDCAFHSGFNEDPFVCEYQQQSSDLPQPPVNAGGGSGGGSGGGSE GGGSEGGGSE GGGSEGGGSGGGSGGDFDYEKMANANKGAMTENADENALQSDAKGKLDS VATDYGAAID GFIGDVSGLANGNGATGDFAGSNSQMAQVGDGDN SPLMNNFRQYLPSLPQS VECRPFVFG AGKPYEFSIDCDKINLFRGVFAFLLYVATFMVVFSTFANILRNKES
Source	in vitro E.coli expression system
Target Names	III
Protein Names	Recommended name: Attachment protein G3P Alternative name(s): Gene 3 protein Short name= G3P Minor coat protein
Expression Region	19-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.