



Recombinant *Oryza sativa* subsp. *japonica* Putative UPF0496 protein 2 (Os06g0718300, LOC_Os06g50410)

Product Code	CSB-CF713692OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5Z8N6
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MIERSNSTPSATPARPPLAVDEEYNQAFRSKSFLLDWSHAHHHLTHTFSSFKL STSTPCA GRGGAREDDFLHAGGDGGAADDSEQSCSYTVLDDFVLEPSPESLARGARLQ QRRRRRPRR HRVETLLIEYFDVTEEACEACSALLAAIGAARRHHLTLRLLLLRLDGGDDDDAK DALARH VRLDNPLSPGSLSEFHDVHARCSPASRLAAAQRRLRLARALRIARGTAAAA LVGACAA AIVAAVLAHALVGIGVAAAAFGATPAGAAARWWGRRAAEKVSSRHYARAGA TLDAAARG AYIVGRDLDTVSRMVRRRAHDELEHGRDVARIAMRGHGERPLLQEVAEEEEEC EEDLRAQL AELEEHVCLCLITINRTRRLVAHEMARGLPPPSPATVTTTSEERLTSS
Source	in vitro E.coli expression system
Target Names	Os06g0718300
Protein Names	Recommended name: Putative UPF0496 protein 2
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.