



Recombinant Chicken Tripartite motif-containing protein 59 (TRIM59)

Product Code	CSB-CF720039CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5ZMD4
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MHQFEEELTCSICYSLFEDPRVLPCSHTEFCRSCLEGVIQLSSNFSIWRPLRVPL KCPNCR SIVEIPASGTESLPINFALKAIIEKYRQEDHSDVATCSEHYRQPLNVYCLLDKKL VCGHC LTIGKHNGHPIDDLHSAYLKEKESSGKILEQLTDKHWSDVCLLIEKLKEQKAQC ESIVQD DKKVVVQYFKKLSETLEHKKQVLLAALDEINRQILEEYEPHIEKLLKIREEQLEL MSLNT SIQKEESPLVFLEKVDNVHQRIKALKEKELPDVKPVEVYPRVGHLLKDVWSKTE IGQINK ILTPKIKLVPKRKLHRSKNSKERGKPEELLQAANPLSVTFIFTVIIAIAVLSFHKPIS SV VIESIPTHISDFFGFLYQDFCTCMQNTVDVVCHKLNSLAEFLGGIVPF
Source	in vitro E.coli expression system
Target Names	TRIM59
Protein Names	Recommended name: Tripartite motif-containing protein 59
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.