



Recombinant Enterobacteria phage I2-2

Attachment protein G3P (III)

Product Code	CSB-CF325353ECP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15415
Product Type	Transmembrane Protein
Immunogen Species	Enterobacteria phage I2-2 (Bacteriophage I2-2)
Sequence	DNWESITKSYYTGFAMSKTVESKDQDGKTVRKEVITQADLTACNDAKASAQD VFNQMKL TFSGIWPDSQFRLVTGDTCVYNGSPSEKTESWSIRAQVEGDMQRSVPDEEPS EQTPEEIC EAKPPIDGVFNNVSKGDEGGFYINYNGCEYEATGVTVCQNDGTVCASSAWKP TGYVPESG ESSSSPVKDGDTGGTGEGGSDTGGDTGGGDTGGGSTGGDTGGSTGGGST GGGSTGGSTGK SLTKEDVTAAIHDA SP SIGDAVKDSL TEDNDQNDNQKKADEQSAKASASVSDA ISDGMRG VGNFVDDLGGESSQYGIGNSEMDLSVSLAKGQLGIDLEGHGS AWESFLNDGA LRPSIPSG HGCTDFVMFQGSVYQLDIGCDKLGDIKSVLSWVMYCLTFWYVFQSATSLLRK GEQ
Source	in vitro E.coli expression system
Target Names	III
Protein Names	Recommended name: Attachment protein G3P Alternative name(s): Gene 3 protein Short name= G3P Minor coat protein
Expression Region	20-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.