



Recombinant *Oryza sativa* subsp. japonica Magnesium transporter MRS2-A, chloroplastic (MRS2-A)

Product Code	CSB-CF859755OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9AUK4
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Sequence	AAGRGGAGLLLLLPPALRAAEGKDGRAVTKDEEEEEAAAAAVEEEGEVEVR REEDKPGD DGSREAAARGSGSGRFSADYISLGIREPVYEVIEVKSNGRMSTKKISRRQLLKS SGLRLR DTRSVDPSLWLMNSMPSLLVREQAILVNLGSLRAIAMHERVLIFNYNSPGGKA FLDSLLP RLNPRNINGGPAMPFQLEVVEAALLSRIQRLERRLMRIEPRVGALLEVLPNRLT ADVLEQ LRLSKQALVELGSRAGDLKQMLIDLLDDPHEIRRICIMGRNCTLDKLSDNMECS VPLEKQ IAEEEEEEIEMLLENYLQRCESIHGQAERLLDSAREMEDSIAVNLSSRRLEVSR VELLLQ VGTFCAIGALIAGIFGMNLKSYLETNAWAFWATTGGIVVGAVAGFFIMYSYLK TRKIL
Source	in vitro E.coli expression system
Target Names	MRS2-A
Protein Names	Recommended name: Magnesium transporter MRS2-A, chloroplastic
Expression Region	56-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.