



# Recombinant Schizosaccharomyces pombe

## Probable glycosidase C21B10.07

### (SPBC21B10.07)

<b>Product Code</b>	CSB-CF891837SXV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9USW3
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Sequence</b>	MGIPDSTTDSRHSLSAALSSASFENIYDPARKNESTNDVIDNHTDTEIDDHDN DHENLD SNNNNENNEAFNEKAAEKLLPWYRRYFIWILIFIVALICSVLIGVLGGVLGHRT AVRDR HPSYKAKTYSLVKEYKGTTFDFGDFMNITDPTHGFVQYLDRNSSAKLGLISAN SSNVIM AADSKHNYSSGRPSIRLQSTQYFEHGLFILDLIHLPGCGTWPAFWTLGDDWP NGGEIDI VEGVNVGTSNQVTLHTGDGCEMEDIKRVMTGTALQTNCWVDAPNSYNAGCG VENPSGPSY GEAFNKNGGGVFVLDWRSEGIRSWFFNRSEIPSDITSGSPQPAKWSEPVADF PDTKCDID KMFSKQKILFDLTFCGDWAGSSVYSSAGCPGSCNDFVGNNPHNFTEAYWNIK SLAVYQY
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	SPBC21B10.07
<b>Protein Names</b>	Recommended name: Probable glycosidase C21B10.07 EC= 3.2.1.-
<b>Expression Region</b>	1-419
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.