



Recombinant Human Mitochondrial chaperone BCS1 (BCS1L)

Product Code	CSB-CF897463HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9Y276
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MPLSDFILALKDNPYFGAGFGLVGVGTALALARKGVQLGLVAFRRHYMITLEVP ARDRSY AWLLSWLTRHSTRTQHLSVETSYLQHESGRISTKFEFVPSPGNHFIWYRGKWI RVERSRE MQMIDLQTGTPWESVTFTALGTDRKVFFNILEEARELALQQEEGKTVMYTAGV SEWRPFG YPRRRRPLNSVVLQQGLADRIVRDVQEFIDNPKWYTDRGIPYRRGYLLYGPPG CGKSSFI TALAGELEHSICLLSLTSSLSDDRLNHLLSVAPQQSLVLEDVDAFLSRDLAV ENPVK YQGLGRLTFSGLLNALDGVASTEARIVFMTTNHVDRLDPALIRPGRVDLKEYV GYCSHWQ LTQMFQRFYPGQAPSLAENFAEHVLRATNQISPAQVQGYFMLYKNPVGAIH NAESLRR
Source	in vitro E.coli expression system
Target Names	BCS1L
Protein Names	Recommended name: Mitochondrial chaperone BCS1 Short name= h-BCS1 Alternative name(s): BCS1-like protein
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.