




Recombinant Epstein-Barr virus Protein BDLF2 (BDLF2)

Product Code	CSB-CF315180EEZ
Abbreviation	Recombinant Epstein-Barr virus BDLF2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C742
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Epstein-Barr virus Protein BDLF2 (BDLF2)
Immunogen Species	Epstein-Barr virus (strain AG876) (HHV-4) (Human herpesvirus 4)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MVDEQVAVEHGTVSHTISREEDGVVHERRVLASGERVEVFYKAPAPRPREGR ASTFHDFTVPAAAAAPGPEPEPEPHPPMPIHANGGGGETKTNTQDQNQNQTTR TRTNAKAEERTAEMDDTMASGGQRGAPISADLLSLSSLTGRMAAMAPSWM KSEVCGERMRFKEDVYDGEAETLAEPPRCFMLS FVFIYYCCYLAFLALLAFGF NPLFLPSFMPVGAKVLRGKGRDFGVPLSYGCPTNPFCKVYTLIPAVVINNVTTY PNNTDSHGHHGGFEAAALHVAALFESGCPNLQAVTNRNRTFNVTRASGRVE RRLVQDMQQRVLASAVVVMHHHCHYETYVFDGVGPEFGTIPTPCFKDVLAFR PSLVTNCTAPLKTSVKGPNWSGAAGGMKRKQCRVDRLTDRSFPAYLEEV MY VMVQ
Research Area	Others
Source	in vitro E.coli expression system
Target Names	BDLF2
Expression Region	1-420aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal10xHis-tagged
Mol. Weight	52.2
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.