



Recombinant Escherichia coli Uncharacterized protein yjiN (yjiN)

Product Code	CSB-CF336567ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P39385
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MNKLIELRRAKRLALSLLLIAAATFVVTFLFPPNFWVSGVKAIAEAMVGALADW FAVVA LFRRVPIPIISRHTAIIPRNKDRIGENLGQFVQEKFLDTQSLVALIRRHEPALLIGN WFS QPENARRVGQHLLQIMSGFLELTDDARIQRLLKRAVHRAIDKVDLSGTSALMLE SMTKND RHQVLLDTLIAQLIALLQRDKSRKFIAQQIVRWLESEHPLKAKILPTEWLGEHSA ELVSD AVNSLLDDISRDRAHQIRHAFDRATFALIDKLKNDPEMAARADAVKSYLKEDEA FNRYLS ELWGDLREWLKVDINSEDSRVKERIARAGQWFGETLIADDALRASLNGHLEQA AHRVAPE FSAFLTRHISDTVKSWDARDMSRQIELNIGKDLQFIRVNGTLVGGCIGLILYLLS QLPAL FPLGNF
Source	in vitro E.coli expression system
Target Names	yjiN
Protein Names	Recommended name: Uncharacterized protein yjiN
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.