



Recombinant Putative neutral sphingomyelinase (T27F6.6)

Product Code	CSB-CF528605CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O45870
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MEPSSEQESQESQKVDMEPVFRRESFAALRRNNVAARQLRVVTLNAWCLPQ PWPIGSTDR VHRLNKIGQYMIDELYDIVGLQELWSYYDFVRLSEQVSSVYPYFHYFHSGFTG SGVCVFS RHPIVSTLTNRYSLNGFAHHIHRGDWFGGKVVLTEIEIDGDLRVNFYTTHLHA EYDREN DLYLPHRTAQAFELAQFVRHTARGADVIVTGDLNMEPCDLGFRLLSHAKLFD AWRMSH EVENEDSEGELLKFRGIAKGGTCDRPDNCYTKRALKNVDDSKRIDYMLFKSGR CNVKLEE CEITLNQIPGEDLNYS DHVGLRARFTIDDRFRHEKSVNTWEPNRPLLIEAIGLVA GGERR ARTDRIFFFILAVICLILILGSLFFEVFPMGFAVLRFALTVVG VFFVWQGLIGLTLE RKA LKA AKQAIQQILNN
Source	in vitro E.coli expression system
Target Names	T27F6.6
Protein Names	Recommended name: Putative neutral sphingomyelinase EC= 3.1.4.12
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.