



Recombinant Mouse Bcl-2-like protein 13 (Bcl2l13)

Product Code	CSB-CF002617MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59017
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MASSTTASLGFHYETKYVVLSYLGLLSQEKKQQGSPPGVQLDVAPQSLNPEVL LKLKSEI EEELKTLDKVESEAFSTSTGFDCHTSPVFSPANPESSIEDCLAHLGERVSQLKE PLQKAL QTILSQPVTYEAYRECTVETAVHASGWNKLLVPLVLLQHLLLELTRRGQEPLR MLLQFGV MYLEEHAAEFIIQQGGWGSVFSLEPEEEEEYPGIIAEDSNDIYILPSDNSGQVSP ESPTV TTSWQSESLPVSLASQSWHTESLPVSLGPESWQQIAMDPPEEVKSLDSSGAG EKSENNSS NSDIVHVEKEEVPEEAFPGAAAPLLTQVPTVEAPEMMRAEKTSPTPSVFVELG EEELEAV TARPEAVERAEGAAQLSEERAGSRKKSHTGEAAAVRGAKSGLPAEGKAVLLF GGAAAVAI LAVAVGVALALRRK
Source	in vitro E.coli expression system
Target Names	Bcl2l13
Protein Names	Recommended name: Bcl-2-like protein 13 Short name= Bcl2-L-13 Alternative name(s): Bcl-rambo Protein Mil1
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.