



Recombinant *Oryza sativa* subsp. *japonica* Putative magnesium transporter MRS2-H (MRS2-H)

Product Code	CSB-CF604684OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10S25
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MALPCAFLSAAAAANATSFSSSPESRRRCRSVHRVPSRPRPPLAPPARVMGKG NSKRKAAN TRLWMRLDRRGCEMILCDKSFVARRSGLPARDLRVLSPLLSRSPSILAREKA MVINLEF VRAIVTADEVLVLEPLAQEVLPFVEKLRKHFPLKSLDVDDVSTHMHTENQDGEL AQDVSC YEVEGANHELPFEFQVLDFALEAVCLSYNSTISDLNRSIAIVLDDLMKSVSTRN LERVWS LKSSLTRLLASVQKVRDEVEHILDDNEAMAHLCTARKTKGQKDLLNTILFPETR LCRTHS SIENSTGIRTCVPSDSDAHILDMLLEAYFKQLDGIRNRIFLVRQYIVDTEDYISIQ L DNK RNELLGLQLTLIIASFIAINTFIAAAFAMNIPHRGYHFVIGVPPFGQFVGATSFLC MSIV ILLFTYAWRNRLCT
Source	in vitro E.coli expression system
Target Names	MRS2-H
Protein Names	Recommended name: Putative magnesium transporter MRS2-H
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.