



Recombinant Uncharacterized protein B0310.3 (B0310.3)

Product Code	CSB-CF608854CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10939
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MPIKCMFPAVNILKFLYHTPRSFLFILDEKKFSNQMEGNFLNNPWINCCFVTT YAGVVA QYFKKGAEFNFNFTTPSIGTFTNEQNEEDSSNYSTSGYDSSAETISANSSPINR SGVRSRI SQKQRQRILKEAHFKAQQLNRKMVVQKSCPPDHEIKPVPSKFYQFDAITDFGF GGPVLLV GMQKDVEVMKLETKEKARSGRKKNRKSKYKCNMYKMTKLAQIVAKIPKKKEVI EIDEDGF QKVSSKKAALRTLKPADVPTPPTKVVENKEEVIKLEVIEQPEPIVLPVSTPTVT FSRFE EMKRVVKVEKAQESAKTKALKKSKAISISRHVGFLLQILEKLEETEEKPQVENEK KVVVKH VQSARKNQKKGRKNRKVEPPKEEFEPYEEDHYNFKRYFLIIGVYVLVFIYVCTN VLTGVGV SYEFPYITLANKSEKNWDSLSKCV
Source	in vitro E.coli expression system
Target Names	B0310.3
Protein Names	Recommended name: Uncharacterized protein B0310.3
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.