



Recombinant Schizosaccharomyces pombe Uncharacterized protein C119.16c (SPBC119.16c)

Product Code	CSB-CF524261SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O42907
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MSIQAIVLATFDAKEGYNVENYYPGDFNVEGIEYLLFPSGIQELDNCTIFFRFQD QLCLS VFSKLQHPSFERSAFFTSVGLILSDDINFGEAVVKYGETLLYIANGLSLATLKYK FGEDA SETYASEKCTSHQLSDSDFFKSLQTSAVNLEFDSLFEKLQGNKFAILGANSKEL SQSYAT ILLDHLGPAFYCLYKFALQRKRILLISSHDDQLYSIIDMIVRLSSIKRSSDASIPILL SD LHPFYSVGLANTSTLLDNDLEEGWIACTTDTVLLSKSSLYDLALYWPDNSFNAN KYPQIF NSNSIRIKPSYDDLINFKGLSRYLSFDGESSWGLTTYSLASKYIFNTSHHNLTDQ EFLNE NMLDYFQRYNQKLLTVLSSNAESFNVSDMQTLGLNPCHSLDKSFVSEISQIWL KKHINWQ YGKYFWLRRVSLIFLASTCFLFILWKLL
Source	in vitro E.coli expression system
Target Names	SPBC119.16c
Protein Names	Recommended name: Uncharacterized protein C119.16c
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.