



Recombinant Escherichia coli Regulator of sigma E protease (rseP)

Product Code	CSB-CF365147ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AEH1
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MLSFLWDLASFIVALGVLITVHEFGHFWARRCGVRVERFSIGFGKALWRRTD KLGTEYV IALIPLGGYVKMLDERAEPVPELRHHAFNNSVQRAAIIAAGPVANFIFAIFAY WLVF IIGVPGVRPVVGEIAANSIAAEAQIAPGTELKAVDGIETPDWDAVRLQLVDKIGD ESTTI TVAPFGSDQRRDVKLDLRHWFEPDKEDPVSSLGIRPRGPQIEPVLENVQPN SAASKAGL QAGDRIVKVDGQPLTQWVTFVMLVRDNP GKSLALEIERQGSPLSLTLIPESKP GNGKAIG FVGIEPKVIPLPDEYKVVRQYGPFN AIVEATDKTWQLMCLTVSMLGKLITGDVK LNNLSG PISIAKGAGMTAELGVVYYLPFLALISVNLGIINLFPLPVL DGGHLLFLAIEKIKGG PVS ERVQDFCYRIGSILLVLLMGLALFNDFSRL
Source	in vitro E.coli expression system
Target Names	rseP
Protein Names	Recommended name: Regulator of sigma E protease EC= 3.4.24.-
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.