



Recombinant *Saccharomyces cerevisiae* Phosphatidylinositol N- acetylglucosaminyltransferase GPI3 subunit (SPT14)

Product Code	CSB-CF461174SVP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3LKQ3
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Sequence	MGFNIAMLCDFYFYPQLGGVEFHLYHLSQKLIDLGHSVVIITHAYKDRVGVRLHTN GLKVY HVPFFVIFRETTTFPTVFSTFPIIRNILLREQIQIVHSHGSASTFAHEGILHANTMGL RTV FTDHSLYGFNNLTSIWWNKLLTFTLTNIDRVICVSNTCKENMIVRTELSPDIISVIP NAV VSEDFKPRDPTGGTKRKQSRDKIVIVVIGRLFPNKGSDLLTRIIPKVCSSHEDVE FIVAG DGPKFIDFQQMIESHRLQKRVLQLLGSPHEKVRDVLCCGDIYLHASLTEAFGTI LVEAAS CNLLIVTTQVGGIPEVLPNEMTVYAEQTSVSDLVQATNKAINIIRSKALDTSSFH DSVSK MYDWMDVAKRTVEIYTNISSTSSADDKDWKMKVANLYKRDGIWAKHLYLLCGI VEYMLFF LLEWLYPRDEIDLAPKWPKKTVSNETKEARET
Source	in vitro E.coli expression system
Target Names	SPT14
Protein Names	Recommended name: Phosphatidylinositol N-acetylglucosaminyltransferase GPI3 subunit EC= 2.4.1.198 Alternative name(s): GlcNAc-PI synthesis protein
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.