



Recombinant *Saccharomyces cerevisiae* Phosphatidylinositol N- acetylglucosaminyltransferase GPI3 subunit (SPT14)

Product Code	CSB-CF335720SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32363
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MGFNIAMLCDFYFYPQLGGVEFHLYHLSQKLIDLGHSVVIITHAYKDRVGVRLHTN GLKVY HVPFFVIFRETTFFPTVFSTFPIIRNILLREQIQIVHSHGSASTFAHEGILHANTMGL RTV FTDHSLYGFNNLTSIWWNKLLTFTLTNIDRVICVSNTCKENMIVRTELSPDIISVIP NAV VSEDFKPRDPTGGTKRKQSRDKIVIVVIGRLFPNKGSDLLTRIIPKVCSSHEDVE FIVAG DGPKFIDFQQMIESHRLQKRQVLLGSPHEKVRDVLCCGDIYLHASLTEAFGTI LVEAAS CNLLIVTTQVGGIPEVLPNEMTVYAEQTSVSDLVQATNKAINIIRSKALDTSSFH DSVSK MYDWMMDVAKRTVEIYTNISSTSSADDKDWKMKVANLYKRDGIWAKHLYLLCGI VEYMLFF LLEWLYPRDEIDLAPKWPKKTVSNETKEARET
Source	in vitro E.coli expression system
Target Names	SPT14
Protein Names	Recommended name: Phosphatidylinositol N-acetylglucosaminyltransferase GPI3 subunit EC= 2.4.1.198 Alternative name(s): GlcNAc-PI synthesis protein
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.