



Recombinant Human Protein odr-4 homolog (ODR4)

Product Code	CSB-CF725554HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5SWX8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGRTYIVEETVGQYLSNINLQGKAFVSGLLIGQCSSQKDYVILATRTPPKKEEQS ENLKHP KAKLDNLDEEWATEHACQVSRMLPGLLVLGVFIITLLELANDFQNALRRLMFA VEKSIN RKRLWNFTEEEVSERVTLHICASTKKIFCRTYDIHDPKSSARPADWKYQSGLS SSWLSLE CTVHINIHIPLSATSVSYLEKNTKNGLTRWAKEIENGVYLINGQVKDEDCDLLE GQKKS SRGNTQATSHSFDVRVLTQLLLNSDHRSTATVQICSGSVNLKGAVKCRAYIHS SKPKVKD AVQAVKRDILNTVADRCEMLFEDLLLNEIPEKKDSEKEFHVLPYRVFVPLPGST VMLCDY KFDDESAAEIRDHFMEMLDHTIQIEDLEIAEETNTACMSSSMNSQASLDNTDDE QPKQPI KTTMLLKIQQNIGVIAAFTVAVLAAGISFHYFSD
Source	in vitro E.coli expression system
Target Names	ODR4
Protein Names	Recommended name: Protein odr-4 homolog Short name= hODR-4 Alternative name(s): LAG1-interacting protein Transactivated by transforming growth factor beta protein 1
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.