



Recombinant Human Transmembrane protein FAM155A (FAM155A)

Product Code	CSB-CF008095HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B1AL88
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MTRGAWMCRQYDDGLKIWLAAPRENEKPFIDSERAQKWRLSLASLLFFTVLLS DHLWFCA EAKLTRARDKEHQQQQRQQQQQQQQQRQRQQQQQQRRQQEPSWPALLAS MGESSPAAQAH RLLSASSSPTLPPSPGDGGGGGGGKGNRGKDDRGKALFLGNSAKPVWRLETC YPQGASSGQ CFTVENADAVCARNWSRGAAGGDGQEVRSKHPTPLWNLSDFYLSFCNSYTL WELFSGLSS PNTLNCSLDVVLKEGGEMTTCRQCVEAYQDYDHHAQEKYEEFESVLHKYLQS EEYSVKSC PEDCKIVYKAWLCSQYFEVTQFNCRKTIPTCKQYCLEVQTRCPFILPDNDEVIYG GLSSFI CTGLYETFLTNDPECCDVRREEKSNNPSKGTVEKSGSCHRTSLTVSSATRL CNSRLKLC VLV LILLHTVLTASAAQNTAGLSFGGINTLEENSTNEE
Source	in vitro E.coli expression system
Target Names	FAM155A
Protein Names	Recommended name: Transmembrane protein FAM155A
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.