



Recombinant Human Metaxin-1 (MTX1)

Product Code	CSB-CF621679HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q13505
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLLGGPPRSRSGTSPKGPWSSTGHVQFGKSPQTWPRRTRPRSPEPAAPS GVRGSTWTRR RDTPRRAGPTALSRVYVGHLMGRRPPSPEARGPVPRSSAASRARRSLASPGI SPGPLTAT IGGAVAGGGPRQGRAEAHKEVFPGQRVGKMAAPMELFCWSGGWGLPSVDL DSLAVLTYAR FTGAPLKVHKISNPWQSPSGTLPALRTSHGEVISVPHKIITHLRKEKYNADYDLS ARQGA DTLAFMSLLEEKLLPVLVHTFWIDTKNYVEVTRKWYAEAMPFPLNFFLPGRMQ RQYMERL QLLTGEHRPEDEEELEKELYREARECLTLLSQRLGSQKFFFGDAPASLDAFVF SYLALLL QAKLPSGKLQVHLRGLHNLCACTHILSLYFPWDGAEVPPQRQTPAGPETEE EPYRRRNQ ILSVLAGLAAMVGYALLSGIVSIQRATPARAPGTRTLGMAEEDDEEE
Source	in vitro E.coli expression system
Target Names	MTX1
Protein Names	Recommended name: Metaxin-1 Alternative name(s): Mitochondrial outer membrane import complex protein 1
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.