



Recombinant Mouse Protein FAM134C (Fam134c)

Product Code	CSB-CF887433MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9CQV4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MEEAEGVAAAPGPASGLAFRGRRAMSGSWERDQQVEAAQRTLVEVLGPYEP LLSRVQAAL VWERPARSALWCLGLNAAFWFFALTSRFRVFLAFSLMIIVCIDQWKNKIWPEI NVPRPD ALDNESWGFVHPRLLSVPELCHHVAEVWVSGTIFIRNLLLFFKKQNPGKFCLLS CGVLTFL AMLGRYIPGLLLSYLMLVIIMMWPLAVYHRLWDRAYVRLKPVLQRLDFSVRGY MMSKQRE RQLRRRALHSERATDSHSDSEEELAAFCPQLDDSTVARELAITDSEHSDAEVS CTENGTF NLSRGQTPLTEGSEDLDGHSDPEESFARDLPDFPSINVDPAGLDDEDDTSIGM PSLMYRS PPGAGDTQVLPASRNEAALPELLSSLPGGSNLTSNLASLVSQGMQLALSEA SQTDPSG PPPRRATRGRFLRAPSSDLDTDAEGDDFELLDQSELNQLDPASSRSH
Source	in vitro E.coli expression system
Target Names	Retreg3
Protein Names	Recommended name: Protein FAM134C
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.