



Recombinant *Oryza sativa* subsp. *indica* Putative magnesium transporter MRS2-G (MRS2-G)

Product Code	CSB-CF383756OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2Z9W7
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MGRRSGGRKLPFFASNASTSSSTKRTRSARRLPSLTRPRASSSPSPASPSPP PPSASHPA PPSPPLAVSPAGAGKVGKKKAGARLWMRLDRWGVSETLHLDKGSIRRAGLP PRDLRILG PVFSDSSSILAREKAMVINLEFIRAIVTADEILLDPLTIDVIPFVEQLTHHLPLKNL VC GNGQPGGDDHGEKHDDSPGDQVPRLNEATGAEHELPPFEFQVLELALETVCS SFDVNVSGL ERRATPVLEELTKNVSTRNLDRVRTLKSDLTRLLAHVQKVRDEIEHLLDDNED MAHLYLT RKQLQNQQVEALISSAASNSIVPGGTSLSRLNNSFRRSVSIATSMHLDNDVEDL EMLLEA YFMQLDGIRNRILSVREYIDDTEDYVNIQLDNQRNELIQLQLTLTIASFGIAVNTFI AGA FAMNIQSKLYSIDDGSFFWPFVGGTSSGCFMICIVLLWYARWKLLGP
Source	in vitro E.coli expression system
Target Names	MRS2-G
Protein Names	Recommended name: Putative magnesium transporter MRS2-G
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.