



Recombinant Uncharacterized glycosyltransferase B0361.8 (B0361.8)

Product Code	CSB-CF346130CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53993
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSDTVISLISHSITTVFYLVPLIALIIPFSLYSGFRRKSKTVAFFHPYCNAGGGGE RVL WAAIRTMQKKFPDHKYFVYSGD TDATKEQILLKARQRFGIELDPSNIQFIYLHW RTLVEA RHYKHCTMLFQALAGLILALEAWFRMVPAVFIDSMGYPLSLPAFRLSGSKVVA YVHYPTI SCDMLDVVESRQETFNNSSTIAQSNVLSWGKLTYYRLFACLYWLAGKAAHVG MVNGSWTQ RHITSIWSRRDVSIVYPPCDVEAFLNIESVAESLLEDTKTVRLLSVGQIRPEKNH KLQLE VLHDVKEPLEKMGYNVELCIAGGCRNEEDQERVKMLKNEAEKLDISEQLIWQL NVPYEDL VVELSKALISIHMTMHNEHFGISVVEAMAASTIILSNDSGGPRMDIVKDYEGHCVG YLSIT KEEYVETILKIVEEGLKKRNDTRKYARKSLTRFGAAAFETHWNKEIEKVL
Source	in vitro E.coli expression system
Target Names	algn-11
Protein Names	Recommended name: Uncharacterized glycosyltransferase B0361.8 EC= 2.4.-.-
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.