



Recombinant Mouse Regulator of microtubule dynamics protein 3 (Fam82a2)

Product Code	CSB-CF663665MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3UJU9
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MSRLGALGGSRAGLGLLLGTAAGLGFLCVLYSQRWKRTQRHGRSHSLPNSLD YAQASERG RQVTQFRAIPGEAGDAAILPSLSQEGQEKVLDRLDFVLTSLMALRREVEELQR SLQGLAG EIVGEVRSHIEENQRVARRRRFPFARERSDSTGSSSVYFTASSGAALTDAESE GGYTTAN AESDYERDSDKESGDAEDEVSCETVRMGRKDSLDDVEAASSPAAAAL EEDD SSGRELVQ LVLLQADELHQGSKQDKREGFQLLLNNKLAYGSRQDFLWRLARAYSDMSDLT EEESGKKS YALNGKEEA EAALKKGDESAACHLWYAVLCGQLAEHEGISKRIQSGFSFKEHV DKAIELQ PEDPRGHFLLGRWCYQVSHLNWLEKKTATALFESPLSATVQDALQSFLKAEEL QPGFSKA GRVYISKCYRELGKNSEARKWMKLAQELPDVTNEDSAFQKDL EEEVILG
Source	in vitro E.coli expression system
Target Names	Rmdn3
Protein Names	Recommended name: Regulator of microtubule dynamics protein 3 Short name= RMD-3 Short name= mRMD-3 Alternative name(s): Protein FAM82A2 Protein FAM82C
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.