



Recombinant Mouse Transmembrane protein FAM155B (Fam155b)

Product Code	CSB-CF008096MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2BDP1
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>MFRGAWMWPGKDAAALTICCCCCCWAPRQSDKPCADSERARWRLSLASL LFFTLLADH LWLCAGARPRARELSSAMRPPWGAGRERQVPPRAVLPPPPSPGEPSSASS GTCGPRYSN LTKAAPAAGSGPVCNGVPEPTGLDAACTKLESLQRLFEPPTPAPPLRPPDSSPS RAPEFPS AKKNLLKGHFRNFTLSFCDTYTVWDLGLMDRPSLDCSLDTLLGDLLAVVAS PGSGTWE ACSNCEAYQRLDRHAQEKYDEFDLVLHKYLQAEYSIRSCTKGCKAVYKAWL CSEYFSV TQQECQNWVPCKQYCLEVQTRCPFILPDNEEMVYGGLPGFICTGLMDTSPKR PETKCCDV QWVSCSESEKKFKDSEPPKTHHQQFHHSYFHYYHHQYHHYHPRHEPPSRVS NKPSLLPVS GGSRLSPSRIRLCVLLVILLHTVVSFSSSQSGGGLGLETLPAL EEGLTQEE</p>
Source	in vitro E.coli expression system
Target Names	Fam155b
Protein Names	Recommended name: Transmembrane protein FAM155B Alternative name(s): Protein TED Transmembrane protein 28
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.