



Recombinant Schizosaccharomyces pombe Uncharacterized protein C1E8.03c (SPBC1E8.03c)

Product Code	CSB-CF528573SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O42968
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MVVIANKGALWAYYCKRLLNSVTYMMYPLIRKRTMKKLLLIVGLLLACSTVMRR IPLFHE SFHLPSLDPRASSTTSQKFQEYRSDFLEKLETAEPEDVIMFTAYGLGVHHTNL FMLACD MAKTSDSQIRFLLLTDGTILPEALYDYNRETVSTCPLSFLSYSTGVERLSKELILK DLLS LQFQQALLAISPSVIVTSEHSPLVMFQAINPYLNNNYTHDTVDNALEENSWIT KLDMQ SLQHFRTPRINVVLIVEDGTYKYLLNLMRDLGRDFKNSEEYPHLFIHLFMSENIP NLSSI RANWPQHRLFINLHFNQKDLNLEIENVWTPNDYTYALVVDLQPDSPQLSSNLIT WLKYKI LLIYYHKSSSTYKNNIAAIVPSFDFSNEEAVILSQTINSNIVLFAPVVFQKFQEYM AVRL LNPNFELPESNGIEFAHEDSVLGHSKPSLTEFHAILGLYSLVISYNHFEGSLSNE YV
Source	in vitro E.coli expression system
Target Names	SPBC1E8.03c
Protein Names	Recommended name: Uncharacterized protein C1E8.03c
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.