



Recombinant Human Protein APCDD1-like (APCDD1L)

Product Code	CSB-CF836717HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NCL9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	EAGGSCLRWEPHCQQPLPDRVPSTAILPPRLNGPWISTGCEVRPGPEFLTRA YTFYPSRL FRAHQFYEDPFCGEPAHSLLVKGVRLRRASWVTRGATEADYHLHKVGIVF HSRRALVD VTGRLNQTRAGRDCARRLPPARAWLPGALYELRSARAQGDCLEALGLTMHEL SLVRVQRR LQPQPRASPRLVEELYLGDHITDPAERRHYRPTGYQRPLQSALHHVQPCPAC GLIARSDV HHPVLPPPLALPLHLGGWWVSSGCEVRPAVLFLTRLFTFHGHSRSWEGYYH HFSDPACR QPTFTVYAAGRYTRGTPSTRVGGTELVEVTRAHVTPMDQVTTAMLNFSEP SSCGGAGA WSMGTERDVTATNGCLPLGIRLPHVEYELFKMEQDPLGQSLLFIGQRPTDGSS PDTPEKR PTSYQAPLVLCHGEAPDFSRPPQHRPSLQKHPSTGGLHIAPFLLPLVLGLAFL HWL
Source	in vitro E.coli expression system
Target Names	APCDD1L
Protein Names	Recommended name: Protein APCDD1-like Alternative name(s): Adenomatosis polyposis coli down-regulated 1 protein-like
Expression Region	25-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.