



Recombinant Rat Sterol regulatory element-binding protein 1 (Srebf1)

Product Code	CSB-CF022657RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56720
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MDELPGEAALQALAEVCEMDAALLTDIEDMLQLINNQSDFPGLFDAPYAG GETGDTG PSSPGASSPESFSSPASLGSSLEAFLGGPKVTPAPLSPPPSAPTAVKMYPSVP PFSPGPG IKEEPVPLTILQPPAPQPSPGTLPPSFPPPPVQLSPAPVLGYSSLPSGFSGLTLP GNTQQ TPSSLPLGSTPGISPTPLHTQVQSSAAQQPPASAAPRMSTVASQIQQVPVVL QPHFIKA DSLLLTAVKTDGTGATMKTAGINTLAPGTAVQAGPLQTLVSGGTILATVPLVVD DKLPIH RLAAGGKALGSAQSRGEKRTAHNAIEKRYRSSINDKIVELKDLVVGTEAKLNKS AVLRKA IDYIRFLQHSNQKQENLTLRSAHKSKSLKDLVSACGSGGGTDVSMEGMKPE VVETLTP PPSDAGSPSQSSPLSLGSRGSSSSGGSDSEPDPSPAFEDNQVKAQRLPSHSRG MLDRSRLAL</p>
Source	in vitro E.coli expression system
Target Names	Srebf1
Protein Names	Recommended name: Sterol regulatory element-binding protein 1 Short name= SREBP-1 Alternative name(s): Adipocyte determination- and differentiation-dependent factor 1 Short name= ADD1 Sterol regulatory element-binding transcripti
Expression Region	1-480
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.