



Recombinant *Saccharomyces cerevisiae* Nuclear rim protein 1 (NUR1)

Product Code	CSB-CF512092SVN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C8Z6L4
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Sequence	<p>MGSNDLINEAYDDSEVVGEEERESKSAWMKRWYQLLTSPDLQLVINEKLEMIN WDAYAKS LAKPLGNFLTILFFIIRLLQDNLIKPNYYKLNKSGAFDLSKSNKLKEFDYLWEIS SSFQ NSNQFYAFQSWYFVTLRFLNNLFRFTIFILLSLNLYVSCKFMFGYFKTYNLFHLK KEFNS PNLTKHNLKDLSKEYYEDIYKQSLWSMLKHFFRGRDDGPHVNQNEVEIFFQL RKWIPTN FMINLFSFSPTAIVFLSFDVSFTSAIAIVFHQYILDYIITKRFQRSVDDDLILSSA AL QEYEDKHIMVRINQCSNIDLSSAMGTRSKTPRIFTTHSLCGEEIREVYNYEKR EFEALP KMTESVPGSRETRIKDYGGISQVSDNQSHPIGFHYSPRMSPYYRDKVLDNNLA QSSNEN LEKGAFLPNQDQNRPSKLSPLRKTPLSARQKRFEGSEFNVLNKNNDINSILR SPKKKKK YHKR</p>
Source	in vitro E.coli expression system
Target Names	NUR1
Protein Names	Recommended name: Nuclear rim protein 1
Expression Region	1-484
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.