



Recombinant Cardiolipin synthase (cls)

Product Code	CSB-CF837963EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8XCC6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MTTVYTLVSWLAILGYWLLIAGVTLRILMKRRVPSAMAWLLIYILPLVGIIAYLA VGE LHLGKRAERARAMWPSTAKWLNDLKACKHIFAEENSSVAAPLFLCERRQGI AGVKGNQ LQLMTESDDVMQALIRDIQLARHNIEMVFYIWQPGGMADQVAESLMAAARRGI HCRLMLD SAGSVAFFHSPWPELMRNAGIEVVEALKVNLMRVFLRRMDLRQHRKMIMIDN YIAYTGSM NMVDPRYFKQDAGVGQWIDLMARMEGPIATAMGIIYSCDWEIETGKRILPPP DVNIMPF EQASGHTIHTIASGPGFPEDLIHQALLTAAYSAREYLMITTPYFVPSDDLHAICT AAQR GVDVSIILPRKNDSMLVGWASRAFFTELLAAGVKVYQFEGGLLHTKSVLVDGE LSLVGTV NLDMRSLWLNFEITLAIDDKGFGADLAAVQDDYISRSRLLDARLWLKRPLWQR VAERLFY FFSPLL
Source	in vitro E.coli expression system
Target Names	clsA
Protein Names	Recommended name: Cardiolipin synthase Short name= CL synthase EC= 2.7.8.-
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.