



Recombinant Arabidopsis thaliana WPP domain-interacting protein 1 (WIP1)

Product Code	CSB-CF812750DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8GXA4
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MDLESESSALESVDDNVLIIQQSASNVCDDGRSLDNGSCSDESVKLLSTSNSV ELGKPMSF DSPGDGGGAYSPVLKGQGLRKWRRIRRD LVKDTSANMENS KALKRGLSGVA HSHGKQMQF QSPEVEQESQGSVGSVNMLKSSGDGFDILGSSGYDSRFVAGVGFSAGMDLEI DDDRSSKS STVARAPKVIRYEKPMISSGQGGNIRVENSKKHARGESVDFEKENSYSSESLS RKQSGRM MDYNGENGETSMRKDDAGGEGGESINTDNRYSDMDPLTEAINGFLALQDAL EKEVQQFQ EIGNEPMPQHHEQVSEANSPHPEIVTLVNNVEQLENMLEETRSMLEVKESHIR DLESTTN QSKHSWGGTEIVVEDIFRQKIEAEIEYLIYSRSIDNLNSQMKLIDEQESLAEET HETLN KLGRVQTKAANFTNRAQDLQNDICIEITGTIKKRACKITSYVLIQLVLLSTVVLLLL SQLP PEPDTVPT
Source	in vitro E.coli expression system
Target Names	WIP1
Protein Names	Recommended name: WPP domain-interacting protein 1
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.