



# Recombinant Arabidopsis thaliana Protein TIC 55, chloroplastic (TIC55)

<b>Product Code</b>	CSB-CF882888DOA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9SK50
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Sequence</b>	<p>           RSAVAGTAVSDQTEGGDVLLNPEEEKRVEVADYDWTEEWYPLYLTKNVPE            DAPLGLTVY            DRQIVLYKDGEGTLRCYEDRCPHRLAKLSEGQLIDGRLECLYHGWQFEGEGK            CVKIPQLP            ASAKIPKAACVKTYEVKDSQGVVWVWMSTKTPPNPEKLPWFENFARPGFFDI            STTHELPY            DHSILLENLMDPAHVPISHDRDFTAKREDAQPLVFEVTERSNRGFAGTWGRE            KEGGKGS            NLLRFDAPCVLQNNREFEGKDGVKNYFSGFLCRPTGQGKSMILVRFVTKR            SPLVSVLP            QFWFHWQNACKVFEQDMGFLSSQNEVLMKEKVPTKDLYLNLKSSDTWVAEYR            KWMDKVGHG            MPYHFGHRTISLPKVPVVEHAPAGLIAALSASYPAKGGIGTMHAPNLANRYFR            HIIHCR            SCSNVIKSFELWKNILSATAVALTALAILVVSQRQWKAVLLGSAALCSAAAYTCLR            AINLN TNNFIRTHRRL         </p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	TIC55
<b>Protein Names</b>	Recommended name: Protein TIC 55, chloroplastic Alternative name(s): Translocon at the inner envelope membrane of chloroplasts 55 Short name= AtTIC55 Translocon at the inner envelope membrane of chloroplasts 55-II
<b>Expression Region</b>	49-539
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.